

Human Chlamydia Pneumoniae IgM (CPNIgM) ELISA Kit | CPNIgM elisa kit

Human Chlamydia Pneumoniae IgM (CPNIgM) ELISA Kit

CATALOG #
AAA204465

REACTIVITY
Human

PRODUCT SYNONYMS

Chlamydia Pneumoniae IgM (CPNIgM); Human Chlamydia Pneumoniae IgM (CPNIgM) ELISA Kit; CPNIgM elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	4ng/ml-64ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CPNIgM elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CPN IGM in the sample, This CPN IGM ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CPN IGM concentration. The concentration of CPN IGM in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.