

Human HAMA ELISA Kit | HAMA elisa kit

Human HAMA ELISA Kit

CATALOG #
AAA204498

REACTIVITY
Human


PRODUCT SYNONYMS

HAMA; Human HAMA ELISA Kit; HAMA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sensitivity	1 ng/mL
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HAMA elisa kit

Intended Uses: This HAMA ELISA kit is intended Laboratory for Research use only and is not for use in test or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of HAMA in the sample, this HAMA ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus HAMA concentration. The concentration of HAMA in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.