

Human Chorionic Gonadotropin ELISA Kit | CGB3 elisa kit

Human Chorionic Gonadotropin Antibody ELISA Kit

CATALOG #
AAA204595

REACTIVITY
Human


PRODUCT SYNONYMS

Chorionic Gonadotropin; Human Chorionic Gonadotropin Antibody ELISA Kit; CGB3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	165
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	25ng/ml-400ng/ml
Sensitivity	1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CGB3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CHORIONIC GONADOTROPIN ANTIBODY in the sample, This CHORIONIC GONADOTROPIN ANTIBODY ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CHORIONIC GONADOTROPIN ANTIBODY concentration. The concentration of CHORIONIC GONADOTROPIN ANTIBODY in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	15426252
NCBI GeneID	1082

NCBI Accession #	NP_149439.1
NCBI GenBank Nucleotide #	NP_149439.1
UniProt Accession #	P0DN86
NCBI Official Full Name	chorionic gonadotropin, beta polypeptide 8
NCBI Official Synonym Full Names	chorionic gonadotropin subunit beta 3
NCBI Official Symbol	CGB3
NCBI Official Synonym Symbols	CGB; CGB5; CGB7; CGB8; hCGB
NCBI Protein Information	choriogonadotropin subunit beta 3
UniProt Protein Name	Choriogonadotropin subunit beta 3
UniProt Gene Name	CGB3
UniProt Synonym Gene Names	CGB; CG-beta

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/chorionic-gonadotropin/204595>