

Human T-Complex Associated Testis Expressed 3 ELISA Kit

Human T-Complex Associated Testis Expressed 3 ELISA Kit

CATALOG #
AAA204709

REACTIVITY
Human

PRODUCT SYNONYMS

T-Complex Associated Testis Expressed 3; Human T-Complex Associated Testis Expressed 3 ELISA Kit; T-Complex Associated Testis Expressed 3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	1.25ng/ml-20ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for T-Complex Associated Testis Expressed 3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of T-COMPLEX ASSOCIATED TESTIS EXPRESSED 3 in the sample, This T-COMPLEX ASSOCIATED TESTIS EXPRESSED 3 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus T-COMPLEX ASSOCIATED TESTIS EXPRESSED 3 concentration. The concentration of T-COMPLEX ASSOCIATED TESTIS EXPRESSED 3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.