

Human HLA Class II Histocompatibility Antigen DQ Chain ELISA Kit

Human HLA Class II Histocompatibility Antigen DQ Chain ELISA Kit

CATALOG #
AAA204890

REACTIVITY
Human


PRODUCT SYNONYMS

HLA Class II Histocompatibility Antigen DQ Chain; Human HLA Class II Histocompatibility Antigen DQ Chain ELISA Kit; HLA Class II Histocompatibility Antigen DQ Chain elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HLA Class II Histocompatibility Antigen DQ Chain elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of HLA CLASS II HISTOCOMPATIBILITY ANTIGEN DQ CHAIN in the sample, This HLA CLASS II HISTOCOMPATIBILITY ANTIGEN DQ CHAIN ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus HLA CLASS II HISTOCOMPATIBILITY ANTIGEN DQ CHAIN concentration. The concentration of HLA CLASS II HISTOCOMPATIBILITY ANTIGEN DQ CHAIN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/hla-class-ii-histocompatibility-antigen-dq-chain/204890>