

Human Protease Activated Receptor 2 ELISA Kit | F2RL1 elisa kit

Human Protease Activated Receptor 2 ELISA Kit

CATALOG #
AAA205092

REACTIVITY
Human

PRODUCT SYNONYMS

Protease Activated Receptor 2; Human Protease Activated Receptor 2 ELISA Kit; F2RL1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	397
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for F2RL1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PROTEASE ACTIVATED RECEPTOR 2 in the sample, This PROTEASE ACTIVATED RECEPTOR 2 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PROTEASE ACTIVATED RECEPTOR 2 concentration. The concentration of PROTEASE ACTIVATED RECEPTOR 2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	34577052
NCBI GeneID	2150

NCBI Accession #	NP_005233.3
NCBI GenBank Nucleotide #	NP_005233.3
UniProt Accession #	P55085
NCBI Official Full Name	proteinase-activated receptor 2 preproprotein
NCBI Official Synonym Full Names	F2R like trypsin receptor 1
NCBI Official Symbol	F2RL1
NCBI Official Synonym Symbols	PAR2; GPR11
NCBI Protein Information	proteinase-activated receptor 2
UniProt Protein Name	Proteinase-activated receptor 2
UniProt Gene Name	F2RL1
UniProt Synonym Gene Names	GPR11; PAR2; PAR-2
UniProt Entry Name	PAR2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/protease-activated-receptor-2/205092>