

Human Peroxisome Proliferators Activator Receptors Alpha ELISA Kit | PPARA elisa kit

Human Peroxisome Proliferators Activator Receptors Alpha ELISA Kit

CATALOG #
AAA205104

REACTIVITY
Human

PRODUCT SYNONYMS

Peroxisome Proliferators Activator Receptors Alpha; Human Peroxisome Proliferators Activator Receptors Alpha ELISA Kit; PPARA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	468
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	50 pg/ml-800 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PPARA elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PEROXISOME PROLIFERATORS ACTIVATOR RECEPTORS alpha in the sample, This PEROXISOME PROLIFERATORS ACTIVATOR RECEPTORS alpha ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PEROXISOME PROLIFERATORS ACTIVATOR RECEPTORS alpha concentration. The concentration of PEROXISOME PROLIFERATORS ACTIVATOR RECEPTORS alpha in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

3041727

NCBI GeneID	5465
UniProt Accession #	Q07869
NCBI Official Full Name	Peroxisome proliferator-activated receptor alpha
NCBI Official Synonym Full Names	peroxisome proliferator activated receptor alpha
NCBI Official Symbol	PPARA
NCBI Official Synonym Symbols	PPAR; NR1C1; hPPAR; PPARalpha
NCBI Protein Information	peroxisome proliferator-activated receptor alpha
UniProt Protein Name	Peroxisome proliferator-activated receptor alpha
UniProt Gene Name	PPARA
UniProt Synonym Gene Names	NR1C1; PPAR; PPAR-alpha
UniProt Entry Name	PPARA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/peroxisome-proliferators-activator-receptors-alpha/205104>