

Human Splicing Factor 3B Subunit 3 ELISA Kit | SF3B3 elisa kit

Human Splicing Factor 3B Subunit 3 ELISA Kit

CATALOG #
AAA205105

REACTIVITY
Human


PRODUCT SYNONYMS

Splicing Factor 3B Subunit 3; Human Splicing Factor 3B Subunit 3 ELISA Kit; SF3B3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1217
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	250 pg/ml-4000 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SF3B3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of SPLICING FACTOR 3B SUBUNIT 3 in the sample, This SPLICING FACTOR 3B SUBUNIT 3 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus SPLICING FACTOR 3B SUBUNIT 3 concentration. The concentration of SPLICING FACTOR 3B SUBUNIT 3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	54112121
NCBI GeneID	23450

NCBI Accession #	NP_036558
NCBI GenBank Nucleotide #	NP_036558.3
UniProt Accession #	Q15393
NCBI Official Full Name	splicing factor 3B subunit 3
NCBI Official Synonym Full Names	splicing factor 3b subunit 3
NCBI Official Symbol	SF3B3
NCBI Official Synonym Symbols	RSE1; SAP130; SF3b130; STAF130
NCBI Protein Information	splicing factor 3B subunit 3
UniProt Protein Name	Splicing factor 3B subunit 3
UniProt Gene Name	SF3B3
UniProt Synonym Gene Names	KIAA0017; SAP130; SF3b130; SAP 130
UniProt Entry Name	SF3B3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/splicing-factor-3b-subunit-3/205105>