

Human Mitogen Activated Protein Kinase 8 ELISA Kit | Mapk8 elisa kit

Human Mitogen Activated Protein Kinase 8 ELISA Kit

CATALOG #
AAA205108

REACTIVITY
Human

PRODUCT SYNONYMS

Mitogen Activated Protein Kinase 8; Human Mitogen Activated Protein Kinase 8 ELISA Kit; Mapk8 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	427
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	62.5 pg/ml-1000 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Mapk8 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of MITOGEN ACTIVATED PROTEIN KINASE 8 in the sample, This MITOGEN ACTIVATED PROTEIN KINASE 8 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus MITOGEN ACTIVATED PROTEIN KINASE 8 concentration. The concentration of MITOGEN ACTIVATED PROTEIN KINASE 8 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	694879445
NCBI GeneID	116554

NCBI GenBank Nucleotide #	NP_446281.2
UniProt Accession #	P49185
NCBI Official Full Name	mitogen-activated protein kinase 8
NCBI Official Synonym Full Names	mitogen-activated protein kinase 8
NCBI Official Symbol	Mapk8
NCBI Official Synonym Symbols	JNK
NCBI Protein Information	mitogen-activated protein kinase 8
UniProt Protein Name	Mitogen-activated protein kinase 8
UniProt Gene Name	Mapk8
UniProt Synonym Gene Names	Jnk1; Prkm8; MAP kinase 8; MAPK 8
UniProt Entry Name	MK08_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/mitogen-activated-protein-kinase-8/205108>