

Human Mothers Against Decapentaplegic Homolog 4 (SMAD4) ELISA Kit | SMAD4 elisa kit

Human Mothers Against Decapentaplegic Homolog 4 (SMAD4) ELISA Kit

CATALOG # AAA205155	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

Mothers Against Decapentaplegic Homolog 4 (SMAD4); Human Mothers Against Decapentaplegic Homolog 4 (SMAD4) ELISA Kit; SMAD4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	552
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	1.25ng/ml-20ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SMAD4 elisa kit	Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of SMAD4 in the sample, This SMAD4 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus SMAD4 concentration. The concentration of SMAD4 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	4885457
NCBI GeneID	4089
NCBI GenBank Nucleotide #	NP_005350.1
UniProt Accession #	Q13485
NCBI Official Full Name	mothers against decapentaplegic homolog 4
NCBI Official Synonym Full Names	SMAD family member 4
NCBI Official Symbol	SMAD4
NCBI Official Synonym Symbols	JIP; DPC4; MADH4; MYHRS

NCBI Protein Information	mothers against decapentaplegic homolog 4
UniProt Protein Name	Mothers against decapentaplegic homolog 4
UniProt Gene Name	SMAD4
UniProt Synonym Gene Names	DPC4; MADH4; MAD homolog 4; Mothers against DPP homolog 4; SMAD 4; Smad4; hSMAD4
UniProt Entry Name	SMAD4_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/mothers-against-decapentaplegic-homolog-4-smad4/205155>