

Human Connective Tissue-Activating Peptide III ELISA Kit

Human Connective Tissue-Activating Peptide III ELISA Kit

CATALOG #
AAA205160

REACTIVITY
Human

PRODUCT SYNONYMS

Connective Tissue-Activating Peptide III; Human Connective Tissue-Activating Peptide III ELISA Kit; Connective Tissue-Activating Peptide III elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	91
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	250ng/ml-4000ng/ml
Sensitivity	10 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Connective Tissue-Activating Peptide III elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration ofCONNECTIVE TISSUE-ACTIVATING PEPTIDE III in the sample, This CONNECTIVE TISSUE-ACTIVATING PEPTIDE III ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versusCONNECTIVE TISSUE-ACTIVATING PEPTIDE III concentration. The concentration ofCONNECTIVE TISSUE-ACTIVATING PEPTIDE III in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	209574
NCBI Official Full Name	connective tissue activating peptide III

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/connective-tissue-activating-peptide-iii/205160>