

Human Ubiquitin Protein Ligase E3 Component N-Recognin 2 ELISA Kit | UBR2 elisa kit

Human Ubiquitin Protein Ligase E3 Component N-Recognin 2 ELISA Kit

CATALOG #
AAA205182

REACTIVITY
Human


PRODUCT SYNONYMS

Ubiquitin Protein Ligase E3 Component N-Recognin 2; Human Ubiquitin Protein Ligase E3 Component N-Recognin 2 ELISA Kit; UBR2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	277
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	1ng/ml-16ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for UBR2 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of UBIQUITIN PROTEIN LIGASE E3 COMPONENT N-RECOGNIN 2 in the sample, This UBIQUITIN PROTEIN LIGASE E3 COMPONENT N-RECOGNIN 2 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus UBIQUITIN PROTEIN LIGASE E3 COMPONENT N-RECOGNIN 2 concentration. The concentration of UBIQUITIN PROTEIN LIGASE E3 COMPONENT N-RECOGNIN 2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

119624502

NCBI GeneID	23304
NCBI Official Full Name	ubiquitin protein ligase E3 component n-recognin 2, isoform CRA_d
NCBI Official Synonym Full Names	ubiquitin protein ligase E3 component n-recognin 2
NCBI Official Symbol	UBR2
NCBI Official Synonym Symbols	C6orf133; bA49A4.1; dJ242G1.1; dJ392M17.3
NCBI Protein Information	E3 ubiquitin-protein ligase UBR2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/ubiquitin-protein-ligase-e3-component-n-recognin-2/205182>