

Human Mixed Lineage Kinase Domain Like Protein ELISA Kit | MLKL elisa kit

Human Mixed Lineage Kinase Domain Like Protein ELISA Kit

CATALOG #
AAA205245

REACTIVITY
Human


PRODUCT SYNONYMS

Mixed Lineage Kinase Domain Like Protein; Human Mixed Lineage Kinase Domain Like Protein ELISA Kit; MLKL elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	471
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	1.25ng/ml-20ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MLKL elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of MIXED LINEAGE KINASE DOMAIN LIKE PROTEIN in the sample, This MIXED LINEAGE KINASE DOMAIN LIKE PROTEIN ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus MIXED LINEAGE KINASE DOMAIN LIKE PROTEIN concentration. The concentration of MIXED LINEAGE KINASE DOMAIN LIKE PROTEIN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	22749323
NCBI GeneID	197259
NCBI GenBank Nucleotide #	NP_689862.1
UniProt Accession #	Q8NB16
NCBI Official Full Name	mixed lineage kinase domain-like protein isoform 1
NCBI Official Synonym Full Names	mixed lineage kinase domain like pseudokinase
NCBI Official Symbol	MLKL

NCBI Official Synonym Symbols	hMLKL
NCBI Protein Information	mixed lineage kinase domain-like protein
UniProt Protein Name	Mixed lineage kinase domain-like protein
UniProt Gene Name	MLKL
UniProt Entry Name	MLKL_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/mixed-lineage-kinase-domain-like-protein/205245>