

Human High Affinity Immunoglobulin Epsilon Receptor Subunit Alpha ELISA Kit | FCER1A elisa kit

Human High Affinity Immunoglobulin Epsilon Receptor Subunit Alpha ELISA Kit

CATALOG #
AAA205247

REACTIVITY
Human


PRODUCT SYNONYMS

High Affinity Immunoglobulin Epsilon Receptor Subunit Alpha; Human High Affinity Immunoglobulin Epsilon Receptor Subunit Alpha ELISA Kit; FCER1A elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	25
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	31.25 pg/ml-500 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FCER1A elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of HIGH AFFINITY IMMUNOGLOBULIN EPSILON RECEPTOR SUBUNIT ALPHA in the sample, This HIGH AFFINITY IMMUNOGLOBULIN EPSILON RECEPTOR SUBUNIT ALPHA ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus HIGH AFFINITY IMMUNOGLOBULIN EPSILON RECEPTOR SUBUNIT ALPHA concentration. The concentration of HIGH AFFINITY IMMUNOGLOBULIN EPSILON RECEPTOR SUBUNIT ALPHA in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	18376638
NCBI GeneID	2205
NCBI Official Full Name	Fc epsilon RI alpha-chain, partial
NCBI Official Synonym Full Names	Fc fragment of IgE receptor Ia
NCBI Official Symbol	FCER1A
NCBI Official Synonym Symbols	FCE1A; FcERI
NCBI Protein Information	high affinity immunoglobulin epsilon receptor subunit alpha

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/high-affinity-immunoglobulin-epsilon-receptor-subunit-alpha/205247>