

Human Lymphoid Enhancer-Binding Factor 1 ELISA Kit | LEF1 elisa kit

Human Lymphoid Enhancer-Binding Factor 1 ELISA Kit

CATALOG #
AAA205312

REACTIVITY
Human


PRODUCT SYNONYMS

Lymphoid Enhancer-Binding Factor 1; Human Lymphoid Enhancer-Binding Factor 1 ELISA Kit; LEF1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	371
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	187.5 pg/ml-3000 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LEF1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Human Lymphoid enhancer-binding factor 1 in the sample, This Human Lymphoid enhancer-binding factor 1 Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus Human Lymphoid enhancer-binding factor 1 concentration. The concentration of Human Lymphoid enhancer-binding factor 1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	195222732
NCBI GeneID	51176

NCBI Accession #	NP_001124185.1
NCBI GenBank Nucleotide #	NP_001124185.1
UniProt Accession #	Q9UJU2
NCBI Official Full Name	lymphoid enhancer-binding factor 1 isoform 2
NCBI Official Synonym Full Names	lymphoid enhancer binding factor 1
NCBI Official Symbol	LEF1
NCBI Official Synonym Symbols	LEF-1; TCF10; TCF7L3; TCF1ALPHA
NCBI Protein Information	lymphoid enhancer-binding factor 1
UniProt Protein Name	Lymphoid enhancer-binding factor 1
UniProt Gene Name	LEF1
UniProt Synonym Gene Names	LEF-1; TCF1-alpha
UniProt Entry Name	LEF1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/lymphoid-enhancer-binding-factor-1/205312>