

Human Phosphorylated C-Jun N-Terminal Kinases/Stress-Activated Protein Kinase ELISA Kit | Prl3b1 elisa kit

Human Phosphorylated C-Jun N-Terminal Kinases/Stress-Activated Protein Kinase ELISA Kit

CATALOG # AAA205332	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

Phosphorylated C-Jun N-Terminal Kinases/Stress-Activated Protein Kinase; Human Phosphorylated C-Jun N-Terminal Kinases/Stress-Activated Protein Kinase ELISA Kit; Prl3b1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	222
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Prl3b1 elisa kit	Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration ofPHOSPHORYLATED C-JUN N-TERMINAL KINASES/STRESS-ACTIVATED PROTEIN KINASEin the sample, This PHOSPHORYLATED C-JUN N-TERMINAL KINASES/STRESS-ACTIVATED PROTEIN KINASEELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versusPHOSPHORYLATED C-JUN N-TERMINAL KINASES/STRESS-ACTIVATED PROTEIN KINASEconcentration. The concentration ofPHOSPHORYLATED C-JUN N-TERMINAL KINASES/STRESS-ACTIVATED PROTEIN KINASEin the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	200405
NCBI GeneID	18776
NCBI Official Full Name	placental lactogen II
NCBI Official Synonym Full Names	prolactin family 3, subfamily b, member 1

NCBI Official Symbol	Prl3b1
NCBI Official Synonym Symbols	PL; PL2; Csh2; Ghd8; PI-2; mplII; mPL-II
NCBI Protein Information	prolactin-3B1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/phosphorylated-c-jun-n-terminal-kinases-stress-activated-protein-kinase/205332>