

Human Vitamin K-Dependent Protein S ELISA Kit | PROS1 elisa kit

Human Vitamin K-Dependent Protein S ELISA Kit

CATALOG #
AAA205499

REACTIVITY
Human

PRODUCT SYNONYMS

Vitamin K-Dependent Protein S; Human Vitamin K-Dependent Protein S ELISA Kit; PROS1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	676
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	6.25ng/ml-100ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PROS1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of VITAMIN K-DEPENDENT PROTEIN S in the sample, This VITAMIN K-DEPENDENT PROTEIN S ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus VITAMIN K-DEPENDENT PROTEIN S concentration. The concentration of VITAMIN K-DEPENDENT PROTEIN S in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	192447438
NCBI GeneID	5627

NCBI Accession #	NP_000304.2
NCBI GenBank Nucleotide #	NP_000304.2
UniProt Accession #	P07225
NCBI Official Full Name	vitamin K-dependent protein S isoform 2 preproprotein
NCBI Official Synonym Full Names	protein S
NCBI Official Symbol	PROS1
NCBI Official Synonym Symbols	PSA; PROS; PS21; PS22; PS23; PS24; PS25; THPH5; THPH6
NCBI Protein Information	vitamin K-dependent protein S
UniProt Protein Name	Vitamin K-dependent protein S
UniProt Gene Name	PROS1
UniProt Synonym Gene Names	PROS
UniProt Entry Name	PROS_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/vitamin-k-dependent-protein-s/205499>