

Human Chloride Channel Accessory 1 ELISA Kit | CLCA1 elisa kit

Human Chloride Channel Accessory 1 ELISA Kit

CATALOG #
AAA205652

REACTIVITY
Human

PRODUCT SYNONYMS

Chloride Channel Accessory 1; Human Chloride Channel Accessory 1 ELISA Kit; CLCA1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	914
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	1.25ng/ml-20ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CLCA1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CHLORIDE CHANNEL ACCESSORY 1 in the sample, This CHLORIDE CHANNEL ACCESSORY 1 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CHLORIDE CHANNEL ACCESSORY 1 concentration. The concentration of CHLORIDE CHANNEL ACCESSORY 1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	110611231
NCBI GeneID	1179

NCBI Accession #	NP_001276.2
NCBI GenBank Nucleotide #	NP_001276.2
UniProt Accession #	A8K7I4
NCBI Official Full Name	calcium-activated chloride channel regulator 1
NCBI Official Synonym Full Names	chloride channel accessory 1
NCBI Official Symbol	CLCA1
NCBI Official Synonym Symbols	CACC; GOB5; CACC1; CLCRG1; CaCC-1; hCLCA1; hCaCC-1
NCBI Protein Information	calcium-activated chloride channel regulator 1
UniProt Protein Name	Calcium-activated chloride channel regulator 1
UniProt Gene Name	CLCA1
UniProt Synonym Gene Names	CACC1; hCLCA1; CaCC-1; hCaCC-1
UniProt Entry Name	CLCA1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/chloride-channel-accessory-1/205652>