

Human Regenerating Islet Derived Protein 3 gamma ELISA Kit | REG3G elisa kit

Human Regenerating Islet Derived Protein 3 gamma ELISA Kit

CATALOG #
AAA205683

REACTIVITY
Human


PRODUCT SYNONYMS

Regenerating Islet Derived Protein 3 gamma; Human Regenerating Islet Derived Protein 3 gamma ELISA Kit; REG3G elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	175

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for REG3G elisa kit

Intended Uses: This REG3R ELISA kit is intended Laboratory for Research use only!. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of REG3R in the sample, this REG3R ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus REG3R concentration. The concentration of REG3R in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	56699413
NCBI GeneID	130120
NCBI Accession #	NP_001008388.1
NCBI GenBank Nucleotide #	NP_001008388.1
UniProt Accession #	Q6UW15
NCBI Official Full Name	regenerating islet-derived protein 3-gamma isoform 1 preproprotein

NCBI Official Synonym Full Names	regenerating family member 3 gamma
NCBI Official Symbol	REG3G
NCBI Official Synonym Symbols	PAP1B; PAPIB; PAP IB; PAP-1B; UNQ429; LPPM429; REG III; REG-III
NCBI Protein Information	regenerating islet-derived protein 3-gamma
UniProt Protein Name	Regenerating islet-derived protein 3-gamma
UniProt Gene Name	REG3G
UniProt Synonym Gene Names	PAP1B; REG-3-gamma; PAP-1B; PAP IB; REG III
UniProt Entry Name	REG3G_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/regenerating-islet-derived-protein-3-gamma/205683>