

Human Mitogen Activated Protein Kinase 6 ELISA Kit | MAPK6 elisa kit

Human Mitogen Activated Protein Kinase 6 ELISA Kit

CATALOG #
AAA205685

REACTIVITY
Human


PRODUCT SYNONYMS

Mitogen Activated Protein Kinase 6; Human Mitogen Activated Protein Kinase 6 ELISA Kit; MAPK6 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	721
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	250 pg/ml-4000 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MAPK6 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of MITOGEN ACTIVATED PROTEIN KINASE 6 in the sample, This MITOGEN ACTIVATED PROTEIN KINASE 6 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus MITOGEN ACTIVATED PROTEIN KINASE 6 concentration. The concentration of MITOGEN ACTIVATED PROTEIN KINASE 6 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4506091
NCBI GeneID	5597

NCBI GenBank Nucleotide #	NP_002739.1
UniProt Accession #	Q16659
NCBI Official Full Name	mitogen-activated protein kinase 6
NCBI Official Synonym Full Names	mitogen-activated protein kinase 6
NCBI Official Symbol	MAPK6
NCBI Official Synonym Symbols	ERK3; PRKM6; p97MAPK; HsT17250
NCBI Protein Information	mitogen-activated protein kinase 6
UniProt Protein Name	Mitogen-activated protein kinase 6
UniProt Gene Name	MAPK6
UniProt Synonym Gene Names	ERK3; PRKM6; MAP kinase 6; MAPK 6; ERK-3; p97-MAPK
UniProt Entry Name	MK06_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/mitogen-activated-protein-kinase-6/205685>