

## Human Voltage Dependent Anion Channel Protein ELISA Kit | VDAC1 elisa kit

Human Voltage Dependent Anion Channel Protein ELISA Kit

CATALOG #  
AAA205709

REACTIVITY  
Human

### PRODUCT SYNONYMS

Voltage Dependent Anion Channel Protein; Human Voltage Dependent Anion Channel Protein ELISA Kit; VDAC1 elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Human                                                             |
| Sequence Length | 283                                                               |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 0.75ng/ml-12ng/ml                                                 |
| Sensitivity     | 0.1 ng/ml                                                         |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for VDAC1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of VOLTAGE DEPENDENT ANION CHANNEL PROTEIN in the sample, This VOLTAGE DEPENDENT ANION CHANNEL PROTEIN ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus VOLTAGE DEPENDENT ANION CHANNEL PROTEIN concentration. The concentration of VOLTAGE DEPENDENT ANION CHANNEL PROTEIN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.