

Human Pyruvate Dehydrogenase E2 ELISA Kit

Human Pyruvate Dehydrogenase E2 ELISA Kit

CATALOG #
AAA205763

REACTIVITY
Human

PRODUCT SYNONYMS

Pyruvate Dehydrogenase E2; Human Pyruvate Dehydrogenase E2 ELISA Kit; Pyruvate Dehydrogenase E2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	543
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	12.5 U/L-200 U/L
Sensitivity	1 U/L

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Pyruvate Dehydrogenase E2 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PYRUVATE DEHYDROGENASE E2 in the sample, This PYRUVATE DEHYDROGENASE E2 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PYRUVATE DEHYDROGENASE E2 concentration. The concentration of PYRUVATE DEHYDROGENASE E2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	39574890
NCBI Official Full Name	pyruvate dehydrogenase E2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/pyruvate-dehydrogenase-e2/205763>