

Human Gamma-Aminobutyric Acid Type A Receptor ELISA Kit

Human Gamma-Aminobutyric Acid Type A Receptor ELISA Kit

CATALOG #
AAA205841

REACTIVITY
Human

PRODUCT SYNONYMS

Gamma-Aminobutyric Acid Type A Receptor; Human Gamma-Aminobutyric Acid Type A Receptor ELISA Kit; Gamma-Aminobutyric Acid Type A Receptor elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	312.5 pg/ml-5000 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Gamma-Aminobutyric Acid Type A Receptor elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of GAMMA-AMINO BUTYRIC ACID TYPE A RECEPTOR in the sample, This GAMMA-AMINO BUTYRIC ACID TYPE A RECEPTOR ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus GAMMA-AMINO BUTYRIC ACID TYPE A RECEPTOR concentration. The concentration of GAMMA-AMINO BUTYRIC ACID TYPE A RECEPTOR in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/gamma-aminobutyric-acid-type-a-receptor/205841>