

Human Proteasome Activator Subunit 1 ELISA Kit | PSME3 elisa kit

Human Proteasome Activator Subunit 1 ELISA Kit

CATALOG #
AAA205910

REACTIVITY
Human

PRODUCT SYNONYMS

Proteasome Activator Subunit 1; Human Proteasome Activator Subunit 1 ELISA Kit; PSME3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	254
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.625ng/ml-10ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PSME3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PROTEASOME ACTIVATOR SUBUNIT 1 in the sample, This PROTEASOME ACTIVATOR SUBUNIT 1 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PROTEASOME ACTIVATOR SUBUNIT 1 concentration. The concentration of PROTEASOME ACTIVATOR SUBUNIT 1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	30410794
NCBI GeneID	10197

NCBI Accession #	NP_005780.2
NCBI GenBank Nucleotide #	NP_005780.2
UniProt Accession #	P61289
NCBI Official Full Name	proteasome activator complex subunit 3 isoform 1
NCBI Official Synonym Full Names	proteasome activator subunit 3
NCBI Official Symbol	PSME3
NCBI Official Synonym Symbols	Ki; PA28G; HEL-S-283; PA28gamma; REG-GAMMA; PA28-gamma
NCBI Protein Information	proteasome activator complex subunit 3
UniProt Protein Name	Proteasome activator complex subunit 3
UniProt Gene Name	PSME3
UniProt Synonym Gene Names	REG-gamma; PA28g; PA28gamma
UniProt Entry Name	PSME3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/proteasome-activator-subunit-1/205910>