

Human Cytoskeleton Associated Protein 4 ELISA Kit | CKAP4 elisa kit

Human Cytoskeleton Associated Protein 4 ELISA Kit

CATALOG #
AAA205925

REACTIVITY
Human

PRODUCT SYNONYMS

Cytoskeleton Associated Protein 4; Human Cytoskeleton Associated Protein 4 ELISA Kit; CKAP4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	602
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	2.5ng/ml-40ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CKAP4 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Human Cytoskeleton associated protein 4 in the sample, This Human Cytoskeleton associated protein 4 Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus Human Cytoskeleton associated protein 4 concentration. The concentration of Human Cytoskeleton associated protein 4 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	19920317
NCBI GeneID	10970

NCBI GenBank Nucleotide #	NP_006816.2
UniProt Accession #	Q07065
NCBI Official Full Name	cytoskeleton-associated protein 4
NCBI Official Synonym Full Names	cytoskeleton associated protein 4
NCBI Official Symbol	CKAP4
NCBI Official Synonym Symbols	p63; CLIMP-63; ERGIC-63
NCBI Protein Information	cytoskeleton-associated protein 4
UniProt Protein Name	Cytoskeleton-associated protein 4
UniProt Gene Name	CKAP4
UniProt Synonym Gene Names	Climp-63; p63
UniProt Entry Name	CKAP4_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cytoskeleton-associated-protein-4/205925>