

Human Myeloid Cell Nuclear Differentiation Antigen ELISA Kit | MNDA elisa kit

Human Myeloid Cell Nuclear Differentiation Antigen ELISA Kit

CATALOG #
AAA205960

REACTIVITY
Human

PRODUCT SYNONYMS

Myeloid Cell Nuclear Differentiation Antigen; Human Myeloid Cell Nuclear Differentiation Antigen ELISA Kit; MNDA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	407
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	125 pg/ml-2000 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MNDA elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of MYELOID CELL NUCLEAR DIFFERENTIATION ANTIGEN in the sample, This MYELOID CELL NUCLEAR DIFFERENTIATION ANTIGEN ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus MYELOID CELL NUCLEAR DIFFERENTIATION ANTIGEN concentration. The concentration of MYELOID CELL NUCLEAR DIFFERENTIATION ANTIGEN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

4505227

NCBI GeneID	4332
NCBI Accession #	NP_002423.1
NCBI GenBank Nucleotide #	NP_002423.1
UniProt Accession #	P41218
NCBI Official Full Name	myeloid cell nuclear differentiation antigen
NCBI Official Synonym Full Names	myeloid cell nuclear differentiation antigen
NCBI Official Symbol	MNDA
NCBI Official Synonym Symbols	PYHIN3
NCBI Protein Information	myeloid cell nuclear differentiation antigen
UniProt Protein Name	Myeloid cell nuclear differentiation antigen
UniProt Gene Name	MNDA
UniProt Entry Name	MNDA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/myeloid-cell-nuclear-differentiation-antigen/205960>