

Human Soluble Urokinase-Type Plasminogen Activator Receptor ELISA Kit | PLAUR elisa kit

Human Soluble Urokinase-Type Plasminogen Activator Receptor ELISA Kit

CATALOG #
AAA206015

REACTIVITY
Human


PRODUCT SYNONYMS

Soluble Urokinase-Type Plasminogen Activator Receptor; Human Soluble Urokinase-Type Plasminogen Activator Receptor ELISA Kit; PLAUR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	301
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	5ng/ml-80ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PLAUR elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of SOLUBLE UROKINASE-TYPE PLASMINOGEN ACTIVATOR RECEPTOR in the sample, This SOLUBLE UROKINASE-TYPE PLASMINOGEN ACTIVATOR RECEPTOR ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus SOLUBLE UROKINASE-TYPE PLASMINOGEN ACTIVATOR RECEPTOR concentration. The concentration of SOLUBLE UROKINASE-TYPE PLASMINOGEN ACTIVATOR RECEPTOR in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	13641309
NCBI GeneID	5329
NCBI Official Full Name	soluble urokinase plasminogen activator receptor
NCBI Official Synonym Full Names	plasminogen activator, urokinase receptor
NCBI Official Symbol	PLAUR
NCBI Official Synonym Symbols	CD87; UPAR; URKR; U-PAR
NCBI Protein Information	urokinase plasminogen activator surface receptor

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/soluble-urokinase-type-plasminogen-activator-receptor/206015>