

Human N-Myc-Interactor (NMI) ELISA Kit | NMI elisa kit

Human N-Myc-Interactor (NMI) ELISA Kit

CATALOG #
AAA206017

REACTIVITY
Human

PRODUCT SYNONYMS

N-Myc-Interactor (NMI); Human N-Myc-Interactor (NMI) ELISA Kit; NMI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	307
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	93.75 pg/ml-1500 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NMI elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of NMI in the sample, This NMIELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus NMI concentration. The concentration of NMI in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	188219599
NCBI GeneID	9111
NCBI GenBank Nucleotide #	NP_004679.2

UniProt Accession #	Q13287
NCBI Official Full Name	N-myc-interactor
NCBI Official Synonym Full Names	N-myc and STAT interactor
NCBI Official Symbol	NMI
NCBI Protein Information	N-myc-interactor
UniProt Protein Name	N-myc-interactor
UniProt Gene Name	NMI
UniProt Synonym Gene Names	Nmi
UniProt Entry Name	NMI_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/n-myc-interactor-nmi/206017>