

Human Casein Kappa ELISA Kit | CSN3 elisa kit

Human Casein Kappa ELISA Kit

CATALOG #
AAA206070

REACTIVITY
Human

PRODUCT SYNONYMS

Casein Kappa; Human Casein Kappa ELISA Kit; CSN3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	180
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	12.5ng/ml-200ng/ml
Sensitivity	1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CSN3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CASEIN KAPPA in the sample, This CASEIN KAPPA ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CASEIN KAPPA concentration. The concentration of CASEIN KAPPA in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	130508685
NCBI Accession #	NP_001076094.1

NCBI GenBank Nucleotide #	NP_001076094.1
UniProt Accession #	P33618
NCBI Official Full Name	kappa-casein
UniProt Protein Name	Kappa-casein
UniProt Gene Name	CSN3
UniProt Synonym Gene Names	CSN10; CSNK
UniProt Entry Name	CASK_RABIT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/casein-kappa/206070>