

Human Proteasome Subunit Alpha Type-3 ELISA Kit | PSMA3 elisa kit

Human Proteasome Subunit Alpha Type-3 ELISA Kit

CATALOG #
AAA206110

REACTIVITY
Human


PRODUCT SYNONYMS

Proteasome Subunit Alpha Type-3; Human Proteasome Subunit Alpha Type-3 ELISA Kit; PSMA3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	248
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.625ng/ml-10ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PSMA3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PROTEASOME SUBUNIT ALPHA TYPE-3 in the sample, This PROTEASOME SUBUNIT ALPHA TYPE-3 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PROTEASOME SUBUNIT ALPHA TYPE-3 concentration. The concentration of PROTEASOME SUBUNIT ALPHA TYPE-3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	13528948
NCBI GeneID	5684

NCBI Official Full Name	Proteasome (prosome, macropain) subunit, alpha type, 3
NCBI Official Synonym Full Names	proteasome 20S subunit alpha 3
NCBI Official Symbol	PSMA3
NCBI Official Synonym Symbols	HC8; PSC3
NCBI Protein Information	proteasome subunit alpha type-3

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/proteasome-subunit-alpha-type-3/206110>