

Human Interferon-Induced 35 KDa Protein ELISA Kit | IFI35 elisa kit

Human Interferon-Induced 35 KDa Protein ELISA Kit

CATALOG #
AAA206166

REACTIVITY
Human

PRODUCT SYNONYMS

Interferon-Induced 35 KDa Protein; Human Interferon-Induced 35 KDa Protein ELISA Kit; IFI35 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	288
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	100 pg/ml-1600 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IFI35 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of INTERFERON-INDUCED 35 KDA PROTEIN in the sample, This INTERFERON-INDUCED 35 KDA PROTEIN ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus INTERFERON-INDUCED 35 KDA PROTEIN concentration. The concentration of INTERFERON-INDUCED 35 KDA PROTEIN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	306922358
NCBI GeneID	3430

NCBI Accession #	NP_005524.2
NCBI GenBank Nucleotide #	NP_005524.2
UniProt Accession #	P80217
NCBI Official Full Name	interferon-induced 35 kDa protein isoform 1
NCBI Official Synonym Full Names	interferon induced protein 35
NCBI Official Symbol	IFI35
NCBI Official Synonym Symbols	IFP35
NCBI Protein Information	interferon-induced 35 kDa protein
UniProt Protein Name	Interferon-induced 35 kDa protein
UniProt Gene Name	IFI35
UniProt Synonym Gene Names	IFP35; IFP 35; Ifi-35
UniProt Entry Name	IN35_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/interferon-induced-35-kda-protein/206166>