

Mouse Peroxisome Proliferators Activator Receptors Gamma (PPAR-Gamma) ELISA Kit | PPAR-gamma elisa kit

Mouse Peroxisome Proliferators Activator Receptors Gamma (PPAR-Gamma) ELISA Kit

CATALOG # AAA206673	REACTIVITY Mouse
-------------------------------	----------------------------


PRODUCT SYNONYMS

Peroxisome Proliferators Activator Receptors Gamma (PPAR-Gamma); Mouse Peroxisome Proliferators Activator Receptors Gamma (PPAR-Gamma) ELISA Kit; PPAR-gamma elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	505
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	75 pg/ml-1200 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PPAR-gamma elisa kit	Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PPAR-gamma in the sample, This PPAR-gamma ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PPAR-gamma concentration. The concentration of PPAR-gamma in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	13432234
NCBI GeneID	5468
UniProt Accession #	P37231
NCBI Official Full Name	Peroxisome proliferator-activated receptor gamma
NCBI Official Synonym Full Names	peroxisome proliferator activated receptor gamma
NCBI Official Symbol	PPARG

NCBI Official Synonym Symbols	GLM1; CIMT1; NR1C3; PPARG1; PPARG2; PPARGgamma
NCBI Protein Information	peroxisome proliferator-activated receptor gamma
UniProt Protein Name	Peroxisome proliferator-activated receptor gamma
UniProt Gene Name	PPARG
UniProt Synonym Gene Names	NR1C3; PPAR-gamma
UniProt Entry Name	PPARG_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/peroxisome-proliferators-activator-receptors-gamma-ppar-gamma/206673>