

Mouse Tissue Inhibitors Of Metalloproteinase 2 (TIMP-2) ELISA Kit | TIMP-2 elisa kit

Mouse Tissue Inhibitors Of Metalloproteinase 2 (TIMP-2) ELISA Kit

CATALOG # AAA206709	REACTIVITY Mouse
-------------------------------	----------------------------


PRODUCT SYNONYMS

Tissue Inhibitors Of Metalloproteinase 2 (TIMP-2); Mouse Tissue Inhibitors Of Metalloproteinase 2 (TIMP-2) ELISA Kit; TIMP-2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	218
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	2.5ng/ml-40ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TIMP-2 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TIMP-2 in the sample, This TIMP-2 Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TIMP-2 concentration. The concentration of TIMP-2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	315585128
NCBI Official Full Name	tissue inhibitor of metalloproteinase 2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/tissue-inhibitors-of-metalloproteinase-2-timp-2/206709>