

## Mouse Gastrointestinal Cancer Marker-CA199 ELISA Kit

Mouse Gastrointestinal Cancer Marker-CA199 ELISA Kit

CATALOG #  
**AAA207243**

REACTIVITY  
**Mouse**

### PRODUCT SYNONYMS

Gastrointestinal Cancer Marker-CA199; Mouse Gastrointestinal Cancer Marker-CA199 ELISA Kit; Gastrointestinal Cancer Marker-CA199 elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Mouse                                                             |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 0.5 U/mL-8 U/mL                                                   |
| Sensitivity     | 0.1 U/ml                                                          |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for Gastrointestinal Cancer Marker-CA199 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of GASTROINTESTINAL CANCER MARKER-CA199 in the sample, This GASTROINTESTINAL CANCER MARKER-CA199 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus GASTROINTESTINAL CANCER MARKER-CA199 concentration. The concentration of GASTROINTESTINAL CANCER MARKER-CA199 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.