

Human Beta-1,4-Galactosyltransferase 1 (b4GALT1) ELISA Kit | b4GALT1 elisa kit

Human Beta-1,4-Galactosyltransferase 1 (b4GALT1) ELISA Kit

CATALOG #
AAA20744

REACTIVITY
Human


PRODUCT SYNONYMS

Beta-1; 4-Galactosyltransferase 1 (b4GALT1); Human Beta-1; 4-Galactosyltransferase 1 (b4GALT1) ELISA Kit; GTB; B4-GALT1; GGTB2; GT1; Beta4Gal-T1; Lactose Synthase A Protein; N-Acetylglucosamine Synthase; Nal synthase; beta-N-acetylglucosamine beta-1; 4-galactosyltransferase 1; b4GALT1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of b4GALT1. No significant cross-reactivity or interference between b4GALT1 and analogues was observed.
Sequence Length	385
Assay Type	Quantitative Sandwich
Samples	Human serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Detection Range	1.56-100ng/mL
Sensitivity	< 0.55ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level b4GALT1 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level b4GALT1 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for b4GALT1 elisa kit	Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to b4GALT1. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to b4GALT1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain b4GALT1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of b4GALT1 in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of b4GALT1 in human serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Product Categories/Family for b4GALT1 elisa kit	Signal transduction; Enzyme & Kinase; Metabolic pathway

NCBI AND UNIPROT INFORMATION

NCBI GI #	13929462
NCBI GeneID	2683
NCBI Accession #	NP_001488.2
NCBI GenBank Nucleotide #	NM_001497.3
UniProt Accession #	P15291; Q12909; Q12910; Q12911; Q14456; Q14509; Q14523; B2R710; D3DRL2
Molecular Weight	42,538 Da
NCBI Official Full Name	beta-1,4-galactosyltransferase 1
NCBI Official Synonym Full Names	UDP-Gal:betaGlcNAc beta 1,4- galactosyltransferase, polypeptide 1
NCBI Official Symbol	B4GALT1
NCBI Official Synonym Symbols	GT1; GTB; CDG2D; GGTB2; B4GAL-T1; beta4Gal-T1
NCBI Protein Information	beta-1,4-galactosyltransferase 1
UniProt Protein Name	Beta-1,4-galactosyltransferase 1
UniProt Gene Name	B4GALT1
UniProt Synonym Gene Names	GGTB2; Beta-1,4-GalTase 1; Beta4Gal-T1; b4Gal-T1
UniProt Entry Name	B4GT1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.