

Mouse Microsomal Triglyceride Transfer Protein Large Subunit ELISA Kit | MTTP elisa kit

Mouse Microsomal Triglyceride Transfer Protein Large Subunit ELISA Kit

CATALOG #
AAA207464

REACTIVITY
Mouse


PRODUCT SYNONYMS

Microsomal Triglyceride Transfer Protein Large Subunit; Mouse Microsomal Triglyceride Transfer Protein Large Subunit ELISA Kit; MTTP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	887
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	93.75 pg/ml-1500 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MTTP elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of MICROSOMAL TRIGLYCERIDE TRANSFER PROTEIN LARGE SUBUNIT in the sample, This MICROSOMAL TRIGLYCERIDE TRANSFER PROTEIN LARGE SUBUNIT ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus MICROSOMAL TRIGLYCERIDE TRANSFER PROTEIN LARGE SUBUNIT concentration. The concentration of MICROSOMAL TRIGLYCERIDE TRANSFER PROTEIN LARGE SUBUNIT in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

1709166

NCBI GeneID	280868
UniProt Accession #	P55156
NCBI Official Full Name	Microsomal triglyceride transfer protein large subunit
NCBI Official Synonym Full Names	microsomal triglyceride transfer protein
NCBI Official Symbol	MTTP
NCBI Official Synonym Symbols	MTP
NCBI Protein Information	microsomal triglyceride transfer protein large subunit
UniProt Protein Name	Microsomal triglyceride transfer protein large subunit
UniProt Gene Name	MTTP
UniProt Synonym Gene Names	MTP
UniProt Entry Name	MTP_BOVIN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/microsomal-triglyceride-transfer-protein-large-subunit/207464>