

Mouse Platelet Endothelial Cell Adhesion Molecule ELISA Kit | Pecam1 elisa kit

Mouse Platelet Endothelial Cell Adhesion Molecule ELISA Kit

CATALOG #
AAA207520

REACTIVITY
Mouse

PRODUCT SYNONYMS

Platelet Endothelial Cell Adhesion Molecule; Mouse Platelet Endothelial Cell Adhesion Molecule ELISA Kit; Pecam1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	678
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	40 pg/ml-640 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Pecam1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PLATELET ENDOTHELIAL CELL ADHESION MOLECULE in the sample, This PLATELET ENDOTHELIAL CELL ADHESION MOLECULE ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PLATELET ENDOTHELIAL CELL ADHESION MOLECULE concentration. The concentration of PLATELET ENDOTHELIAL CELL ADHESION MOLECULE in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	84029389
NCBI GeneID	29583

UniProt Accession #	Q3SWT0
NCBI Official Full Name	Platelet endothelial cell adhesion molecule
NCBI Official Synonym Full Names	platelet and endothelial cell adhesion molecule 1
NCBI Official Symbol	Pecam1
NCBI Official Synonym Symbols	CD31; Pecam
NCBI Protein Information	platelet endothelial cell adhesion molecule
UniProt Protein Name	Platelet endothelial cell adhesion molecule
UniProt Gene Name	Pecam1
UniProt Synonym Gene Names	Pecam; PECAM-1
UniProt Entry Name	PECA1_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/platelet-endothelial-cell-adhesion-molecule/207520>