

Mouse Urinary Type Plasminogen Activator ELISA Kit

Mouse Urinary Type Plasminogen Activator ELISA Kit

CATALOG #
AAA207521

REACTIVITY
Mouse

PRODUCT SYNONYMS

Urinary Type Plasminogen Activator; Mouse Urinary Type Plasminogen Activator ELISA Kit; Urinary Type Plasminogen Activator elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	6.25ng/ml-100ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Urinary Type Plasminogen Activator elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of URINARY TYPE PLASMINOGEN ACTIVATOR in the sample, This URINARY TYPE PLASMINOGEN ACTIVATOR ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus URINARY TYPE PLASMINOGEN ACTIVATOR concentration. The concentration of URINARY TYPE PLASMINOGEN ACTIVATOR in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.