

Mouse Phosphoinositide-Dependent Protein Kinase-1 ELISA Kit

Mouse Phosphoinositide-Dependent Protein Kinase-1 ELISA Kit

CATALOG #
AAA207566

REACTIVITY
Mouse


PRODUCT SYNONYMS

Phosphoinositide-Dependent Protein Kinase-1; Mouse Phosphoinositide-Dependent Protein Kinase-1 ELISA Kit; Phosphoinositide-Dependent Protein Kinase-1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	564
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.625ng/ml-10ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Phosphoinositide-Dependent Protein Kinase-1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PHOSPHOINOSITIDE-DEPENDENT PROTEIN KINASE-1 in the sample, This PHOSPHOINOSITIDE-DEPENDENT PROTEIN KINASE-1 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PHOSPHOINOSITIDE-DEPENDENT PROTEIN KINASE-1 concentration. The concentration of PHOSPHOINOSITIDE-DEPENDENT PROTEIN KINASE-1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

345050324

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/phosphoinositide-dependent-protein-kinase-1/207566>