

Mouse Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit | UCHL1 elisa kit

Mouse Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit

CATALOG #
AAA20777

REACTIVITY
Mouse


PRODUCT SYNONYMS

Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1); Mouse Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit; UCHL1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1). No significant cross-reactivity or interference between Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) and analogues was observed.
Sequence Length	223
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, Cell lysates, Cell culture supernates and other biological fluids
Detection Range	0.78-50ng/mL
Sensitivity	0.32ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for UCHL1 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Product Categories/Family for UCHL1 elisa kit	Enzyme & Kinase

NCBI AND UNIPROT INFORMATION

NCBI GI #	21361091
NCBI GeneID	7345
NCBI Accession #	NP_004172.2
NCBI GenBank Nucleotide #	NM_004181.4
UniProt Accession #	P09936; Q4W5K6; Q71UM0
Molecular Weight	24,824 Da
NCBI Official Full Name	ubiquitin carboxyl-terminal hydrolase isozyme L1
NCBI Official Synonym Full Names	ubiquitin carboxyl-terminal esterase L1 (ubiquitin thiolesterase)
NCBI Official Symbol	UCHL1
NCBI Official Synonym Symbols	NDGOA; PARK5; PGP95; PGP9.5; Uch-L1; HEL-117; PGP 9.5
NCBI Protein Information	ubiquitin carboxyl-terminal hydrolase isozyme L1; epididymis luminal protein 117; neuron cytoplasmic protein 9.5; ubiquitin C-terminal hydrolase; ubiquitin thioesterase L1
UniProt Protein Name	Ubiquitin carboxyl-terminal hydrolase isozyme L1
UniProt Gene Name	UCHL1
UniProt Synonym Gene Names	UCH-L1; PGP9.5
UniProt Entry Name	UCHL1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

