

Mouse Nicotinamide Adenine Dinucleotide Phosphate Oxidase ELISA Kit

Mouse Nicotinamide Adenine Dinucleotide Phosphate Oxidase ELISA Kit

CATALOG #
AAA207916

REACTIVITY
Mouse


PRODUCT SYNONYMS

Nicotinamide Adenine Dinucleotide Phosphate Oxidase; Mouse Nicotinamide Adenine Dinucleotide Phosphate Oxidase ELISA Kit; Nicotinamide Adenine Dinucleotide Phosphate Oxidase elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	7.5 U/L-120 U/L
Sensitivity	0.1 U/L

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Nicotinamide Adenine Dinucleotide Phosphate Oxidase elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of NICOTINAMIDE ADENINE DINUCLEOTIDE PHOSPHATE OXIDASE in the sample, This NICOTINAMIDE ADENINE DINUCLEOTIDE PHOSPHATE OXIDASE ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus NICOTINAMIDE ADENINE DINUCLEOTIDE PHOSPHATE OXIDASE concentration. The concentration of NICOTINAMIDE ADENINE DINUCLEOTIDE PHOSPHATE OXIDASE in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.