

Mouse Neutrophil Gelatinase-Associated Lipocalin ELISA Kit | LCN2 elisa kit

Mouse Neutrophil Gelatinase-Associated Lipocalin ELISA Kit

CATALOG #
AAA207977

REACTIVITY
Mouse


PRODUCT SYNONYMS

Neutrophil Gelatinase-Associated Lipocalin; Mouse Neutrophil Gelatinase-Associated Lipocalin ELISA Kit; LCN2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	198
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	10 mU/mL-160 mU/mL
Sensitivity	1 mU/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LCN2 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of NEUTROPHIL GELATINASE-ASSOCIATED LIPOCALIN in the sample, This NEUTROPHIL GELATINASE-ASSOCIATED LIPOCALIN ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus NEUTROPHIL GELATINASE-ASSOCIATED LIPOCALIN concentration. The concentration of NEUTROPHIL GELATINASE-ASSOCIATED LIPOCALIN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	1171700
NCBI GeneID	3934
UniProt Accession #	P80188
NCBI Official Full Name	Neutrophil gelatinase-associated lipocalin
NCBI Official Synonym Full Names	lipocalin 2
NCBI Official Symbol	LCN2
NCBI Official Synonym Symbols	p25; 24p3; MSFI; NGAL

NCBI Protein Information	neutrophil gelatinase-associated lipocalin
UniProt Protein Name	Neutrophil gelatinase-associated lipocalin
UniProt Gene Name	LCN2
UniProt Synonym Gene Names	HNL; NGAL; NGAL
UniProt Entry Name	NGAL_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/neutrophil-gelatinase-associated-lipocalin/207977>