

Mouse Histone Deacetylase 6 ELISA Kit | Hdac6 elisa kit

Mouse Histone Deacetylase 6 ELISA Kit

CATALOG #
AAA208189

REACTIVITY
Mouse

PRODUCT SYNONYMS

Histone Deacetylase 6; Mouse Histone Deacetylase 6 ELISA Kit; Hdac6 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	1149
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	12.5 U/L-200 U/L
Sensitivity	1 U/L

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Hdac6 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of HISTONE DEACETYLASE 6 in the sample, This HISTONE DEACETYLASE 6 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus HISTONE DEACETYLASE 6 concentration. The concentration of HISTONE DEACETYLASE 6 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	341941079
NCBI GeneID	15185

UniProt Accession #	Q9Z2V5
NCBI Official Full Name	Histone deacetylase 6
NCBI Official Synonym Full Names	histone deacetylase 6
NCBI Official Symbol	Hdac6
NCBI Official Synonym Symbols	Hd6; Sfc6; Hdac5; mHDA2
NCBI Protein Information	histone deacetylase 6
UniProt Protein Name	Histone deacetylase 6
UniProt Gene Name	Hdac6
UniProt Synonym Gene Names	HD6
UniProt Entry Name	HDAC6_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/histone-deacetylase-6/208189>