

Mouse Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 ELISA Kit

Mouse Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 ELISA Kit

CATALOG #
AAA208226

REACTIVITY
Mouse


PRODUCT SYNONYMS

Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2; Mouse Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 ELISA Kit; Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	125 pg/ml-2000 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TYROSINE KINASE WITH IMMUNOGLOBULIN-LIKE AND EGF-LIKE DOMAINS 2 in the sample, This TYROSINE KINASE WITH IMMUNOGLOBULIN-LIKE AND EGF-LIKE DOMAINS 2 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TYROSINE KINASE WITH IMMUNOGLOBULIN-LIKE AND EGF-LIKE DOMAINS 2 concentration. The concentration of TYROSINE KINASE WITH IMMUNOGLOBULIN-LIKE AND EGF-LIKE DOMAINS 2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.