

## Mouse Lipase, Monoacylglycerol (MGL) ELISA Kit | MGL elisa kit

Mouse Lipase, Monoacylglycerol (MGL) ELISA Kit

CATALOG #  
**AAA208353**

REACTIVITY  
**Mouse**

### PRODUCT SYNONYMS

Lipase; Monoacylglycerol (MGL); Mouse Lipase; Monoacylglycerol (MGL) ELISA Kit; MGL elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Mouse                                                             |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 7.5 U/L-120 U/L                                                   |
| Sensitivity     | 0.1 U/L                                                           |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for MGL elisa kit

**Intended Use:** The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of MGL in the sample, This MGL ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus MGL concentration. The concentration of MGL in the samples is then determined by comparing the O.D. of the samples to the standard curve.

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.