

Mouse Urokinase Type Plasminogen Activator Receptor (UPAR) ELISA Kit | UPAR elisa kit

Mouse Urokinase Type Plasminogen Activator Receptor (UPAR) ELISA Kit

CATALOG #
AAA208362

REACTIVITY
Mouse

PRODUCT SYNONYMS

Urokinase Type Plasminogen Activator Receptor (UPAR); Mouse Urokinase Type Plasminogen Activator Receptor (UPAR) ELISA Kit; UPAR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	62.5 pg/ml-1000 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for UPAR elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of UPAR in the sample, This UPAR ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus UPAR concentration. The concentration of UPAR in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/urokinase-type-plasminogen-activator-receptor-upar/208362>