

Rat Bone Morphogenetic Protein Receptor II (BMPR-II) ELISA Kit | BMPR-II elisa kit

Rat Bone Morphogenetic Protein Receptor II (BMPR-II) ELISA Kit

CATALOG #
AAA208835

REACTIVITY
Rat

PRODUCT SYNONYMS

Bone Morphogenetic Protein Receptor II (BMPR-II); Rat Bone Morphogenetic Protein Receptor II (BMPR-II) ELISA Kit; BMPR-II elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Sequence Length	1038
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BMPR-II elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of BMPR- II in the sample, This BMPR- II ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus BMPR- II concentration. The concentration of BMPR- II in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	15451916
NCBI GeneID	659

NCBI GenBank Nucleotide #	NP_001195.2
UniProt Accession #	Q13873
NCBI Official Full Name	bone morphogenetic protein receptor type-2
NCBI Official Synonym Full Names	bone morphogenetic protein receptor type 2
NCBI Official Symbol	BMPR2
NCBI Official Synonym Symbols	BMR2; PPH1; BMPR3; BRK-3; POVD1; T-ALK; BMPR-II
NCBI Protein Information	bone morphogenetic protein receptor type-2
UniProt Protein Name	Bone morphogenetic protein receptor type-2
UniProt Gene Name	BMPR2
UniProt Synonym Gene Names	PPH1; BMP type-2 receptor; BMPR-2; BMP type II receptor; BMPR-II
UniProt Entry Name	BMPR2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/bone-morphogenetic-protein-receptor-ii-bmpr-ii/208835>