

Rat Tissue Inhibitors Of Metalloproteinase 4 (TIMP-4) ELISA Kit | TIMP-4 elisa kit

Rat Tissue Inhibitors Of Metalloproteinase 4 (TIMP-4) ELISA Kit

CATALOG #
AAA208902

REACTIVITY
Rat


PRODUCT SYNONYMS

Tissue Inhibitors Of Metalloproteinase 4 (TIMP-4); Rat Tissue Inhibitors Of Metalloproteinase 4 (TIMP-4) ELISA Kit; TIMP-4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Sequence Length	51
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	62.5 pg/ml-1000 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TIMP-4 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TIMP-4 in the sample, This TIMP-4 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TIMP-4 concentration. The concentration of TIMP-4 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	48734400
NCBI GeneID	110595

NCBI Official Full Name	tissue inhibitor of metalloproteinase 4, partial
NCBI Official Synonym Full Names	tissue inhibitor of metalloproteinase 4
NCBI Official Symbol	Timp4
NCBI Official Synonym Symbols	Timp-4
NCBI Protein Information	metalloproteinase inhibitor 4

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/tissue-inhibitors-of-metalloproteinase-4-timp-4/208902>