

Rat Urokinase-Type Plasminogen Activator (uPA/PLAU) ELISA Kit | uPA/PLAU elisa kit

Rat Urokinase-Type Plasminogen Activator (uPA/PLAU) ELISA Kit

CATALOG #
AAA209066

REACTIVITY
Rat


PRODUCT SYNONYMS

Urokinase-Type Plasminogen Activator (uPA/PLAU); Rat Urokinase-Type Plasminogen Activator (uPA/PLAU) ELISA Kit; uPA/PLAU elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Sequence Length	345
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	75 pg/ml-1200 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for uPA/PLAU elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of UPA;PLAU in the sample, This UPA;PLAU ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus UPA;PLAU concentration. The concentration of UPA;PLAU in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	984290529
NCBI GeneID	5328

NCBI GenBank Nucleotide #	NP_001306120.1
UniProt Accession #	P00749
NCBI Official Full Name	urokinase-type plasminogen activator isoform 3
NCBI Official Synonym Full Names	plasminogen activator, urokinase
NCBI Official Symbol	PLAU
NCBI Official Synonym Symbols	ATF; QPD; UPA; URK; u-PA; BDPLT5
NCBI Protein Information	urokinase-type plasminogen activator
UniProt Protein Name	Urokinase-type plasminogen activator
UniProt Gene Name	PLAU
UniProt Synonym Gene Names	U-plasminogen activator; uPA
UniProt Entry Name	UROK_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/urokinase-type-plasminogen-activator-upa-plau/209066>