

## Rat Ketosteroid Isomerase ELISA Kit

Rat Ketosteroid Isomerase ELISA Kit

CATALOG #  
**AAA210155**

REACTIVITY  
**Rat**


### PRODUCT SYNONYMS

Ketosteroid Isomerase; Rat Ketosteroid Isomerase ELISA Kit; Ketosteroid Isomerase elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Rat                                                               |
| Sequence Length | 125                                                               |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 2ng/ml-32ng/ml                                                    |
| Sensitivity     | 0.1 ng/ml                                                         |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for Ketosteroid Isomerase elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Rat ketosteroid isomerase in the sample, This Rat ketosteroid isomerase Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus Rat ketosteroid isomerase concentration. The concentration of Rat ketosteroid isomerase in the samples is then determined by comparing the O.D. of the samples to the standard curve.

### NCBI AND UNIPROT INFORMATION

|                         |                       |
|-------------------------|-----------------------|
| NCBI GI #               | 151322                |
| NCBI Official Full Name | ketosteroid isomerase |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/ketosteroid-isomerase/210155>